



CUS-3000

Digital PalmTop Ultrasound Scanner

Features:

- Advanced digital imaging technologies
- Incredible clarity
- Abundant software packages
- Ergonomic design, user friendly
- Cost effective

CUS-3000

Digital PalmTop Ultrasound Scanner

Specifications:

Scanning:	electronic convex & linear
Mode:	B, B+B, B+M, M, 4B
Gray scale:	256
Scanning depth:	240mm(max)
Cine-loop:	819 frames(max)
Image storage:	permanent 32 frames
Display:	6.4 " LCD
Body mark:	40
Notation:	ID, age, sex, comment, date & time, hospital, doctor, language, OB table
Application:	abdominal, OB, GYN, urology, cardiology, pediatric, small parts, etc.
Measurement:	distance, area, girth, volume, EF rate, HR
OB measurement:	EDD, GA, FW (by measurement of BPD, GS, CRL, FL, HC)
Image conversion:	up/down, right/left, black/white
Local zoom:	x2 in real-time or frozen
Video out:	PAL-D, video printer interface, USB interface
Power supply:	110/220V, 50/60Hz
Battery:	built-in rechargeable, 3 hours' working time
Net weight:	1.3kgs
Gross weight:	6.5kgs



Standard Configurations

CUS-3000 main unit
R40/3.5MHz convex probe
rechargeable battery
AC adaptor
charging cable
metal case
USB drive disk and connecting cable

Option Configurations

7.5MHz endorectal
7.5MHz linear
R20/5.0MHz micro-convex
video printer
car charger



7.5MHz Linear



20R/5.0MHz Micro-Convex



13R/6.5MHz Transvaginal Convex



7.5MHz Endorectal Linear



Carewell Co., Ltd.

Add: 5A, Huating 3, Time City, Chuangye Rd., Nanshan, Shenzhen 518054, China
Tel: +86-755-8617 0472 8617 0684 Fax: +86-755-8617 0478
E-mail: market@carewell.com.cn Website: www.carewell.com.cn



CUS-3000V

Digital PalmTop Ultrasound Scanner for Veterinary

Specifications:

Scanning:	electronic convex & linear
Mode:	B, B+B, B+M, M, 4B
Gray scale:	256
Scanning depth:	60-130mm(rectal probe)
Cine-loop:	508 frames(max)
Image storage:	permanent 32 frames
Display:	6.4 " LCD
Notation:	ID, age, sex, comment, date & time, hospital, doctor
Application:	abdominal, OB, GYN, urology, cardiology, pediatric, small parts, etc.
Measurement:	distance, area, girth, volume, EF rate, HR
OB measurement:	Equine-GA(GSD) Bovine-GA(BL, SL, HL) Swine-GA(HL) Sheep-GA(USD) Cat-GA, EDD(HD, BD) Dog-GA, EDD(GSD, CRL, HD, BD)
Image conversion:	up/down, right/left, black/white
Local zoom:	x2 in real-time or frozen
Video out:	PAL-D, video printer interface, USB interface
Power supply:	110/220V, 50/60Hz
Battery:	built-in rechargeable, 3 hours' working time
Net weight:	1.3kgs
Gross weight:	6.5kgs



Standard Configurations

- CUS-3000V main unit
- 6.5MHz rectal Probe
- Rechargeable battery
- AC adaptor
- Charging cable
- Metal case
- USB drive disk and connecting cable



Option Configurations

- R40MHz/3.5MHz convex
- 7.5MHz endorectal
- 7.5MHz linear
- R20/5.0MHz micro-convex
- Video printer
- Var charger



40R/3.5MHz Convex



7.5MHz Linear



20R/5.0MHz Micro-Convex



7.5MHz Endorectal Linear



Carewell Co., Ltd.

Add: 5A, Huating 3, Time City, Chuangye Rd., Nanshan, Shenzhen 518054, China
Tel: +86-755-8617 0472 8617 0684 Fax: +86-755-8617 0478
E-mail: market@carewell.com.cn Website: www.carewell.com.cn

